

21-MINUTES SESSION

... FOR PREVENTION, REGENERATION,
ANTI-AGING, WELL-BEING,
SLEEP QUALITY AND RELAXATION.

The Airnergy Professional Plus is equipped with 4 Airnergy activation units and colour light therapy (blue). The recommended application time is 21 minutes.



EQUIPMENT

- 4 Airnergy catalysts (activation units)
- Integrated colour light therapy (blue)
- 3 intensity levels (50 %, 75 %, 100 %)
- LCD display
- Five-language menu guidance (DE, EN, FR, IT, ES)
- Microcontroller-controlled electronics
- Individual programme settings
- Automatic self-test
- Switchable acoustic signal
- Optional: aromatherapy

TECHNICAL SPECIFICATIONS

Power supply	12 V DC
Power supply unit	100-240V AC/50-60Hz out: 2.5 A
Power consumption	max. 1,600 mA
Air supply	approx. 4 l per min.
Dimensions (mm) incl. glass bottle	350 x 330 x 350 (L x W x D)
Weight	4.4 kg

INFORMATION

Service	Every 2,000 hours of operation (approx. 6,000 sessions)
Warranty	2 years


Shipment includes	1 power supply unit, 20 nasal cannulas, 1 glass bottle incl. sparkling element, 2 air filters, 1 cannula bag, 1 user manual
-------------------	---

Accessories (optional)	Aroma bottle, stainless steel support for aroma bottle
------------------------	--

Article-No.	2040055
-------------	---------

EAN-No.	0605930524147
---------	---------------





21-MINUTES SESSION

... FOR PREVENTION, REGENERATION,
ANTI-AGING, WELL-BEING,
SLEEP QUALITY AND RELAXATION.

The Airnergy Professional Plus is equipped with 4 Airnergy activation units and colour light therapy (blue). The recommended application time is 21 minutes.

EQUIPMENT

- 4 Airnergy catalysts (activation units)
- Integrated colour light therapy (blue)
- 3 intensity levels (50 %, 75 %, 100 %)
- LCD display
- Five-language menu guidance (DE, EN, FR, IT, ES)
- Microcontroller-controlled electronics
- Individual programme settings
- Automatic self-test
- Switchable acoustic signal
- Optional: aromatherapy

TECHNICAL SPECIFICATIONS

Power supply	12 V DC
Power supply unit	100-240V AC/50-60Hz out: 2.5 A
Power consumption	max. 1,600 mA
Air supply	approx. 4 l per min.
Dimensions (mm) incl. glass bottle	350 x 310 x 240 (L x W x D)
Weight	4.4 kg

INFORMATION

Service	Every 2,000 hours of operation (approx. 6,000 sessions)
Warranty	2 years

Shipment includes	1 power supply unit, 20 nasal cannulas, 1 glass bottle incl. sparkling element, 2 air filters, 1 cannula bag, 1 user manual
Accessories (optional)	Aroma bottle, stainless steel support for aroma bottle
Article-No.	2040055 (compact)
EAN-No.	0605930524147

